


Difração de Raio X

J. Kogler – EP USP

2016



a Campo de conexão à rede

b Campo de comando

c Campo de conexão

d Câmara do tubo
(com tubo de raios X)

e Câmara de experiências
(aqui com goniômetro)


f Tela fluorescente



g Canal livre


h Tecla de bloqueio


i Pés



k Alças de transporte















Bragg Planck Transmission Moseley

$\beta / ^\circ$

R
1/s



X-Ray Apparatus


U = 5.0 kV
I = 0.00 mA
 $\Delta t = 1$ s
 $\Delta \beta = 0.1$ °

0 /s
0.0 °
0.0 °

Salvar / apagar area de trabalho

Parar / reiniciar medições

Ajustes dos parâmetros experimentais



Visualização dos ângulos

2ª parte – determinação de h via Duane-Hunt

1ª parte – difratometria de Bragg p/medir " d "

X-Ray Apparatus




Bragg Planck Transmission Moseley


$\beta / ^\circ$

R 10
1/s


X-Ray Apparatus

U = 5.0 kV
I = 0.00 mA
 $\Delta t = 1$ s
 $\Delta\beta = 0.1^\circ$





Bragg Planck Transmission Moseley



X-Ray Apparatus

U = 5.0 kV
I = 0.00 mA
 $\Delta t = 1$ s
 $\Delta \beta = 0.1^\circ$

Settings

X-Ray Apparatus Crystal General

Goniometer
 Sensor
 Target
 Coupled $\beta = 0.0^\circ$

Parameter
U = 5.0 kV $\beta_{min} = 2.5^\circ$ $\Delta \beta = 0.1^\circ$
I = 0.00 mA $\beta_{max} = 30.0^\circ$ $\Delta t = 1$ s

Current actions
 Scan HV Replay Zero Loudspeaker


Close Crystal calibration Help

X-Ray Apparatus


Bragg Planck Transmission Moseley

$\beta / ^\circ$ R 10
1/s



Bragg Planck Transmission Moseley



Bragg Planck Transmission Moseley



Bragg Planck Transmission Moseley



Settings

X-Ray Apparatus Crystal General

Goniometer

Sensor
 Target $\beta = 6.9^\circ$
 Coupled

Parameter

U = 35.0 kV $\beta_{\min} = 2.5^\circ$ $\Delta\beta = 0.1^\circ$
I = 1.00 mA $\beta_{\max} = 30.0^\circ$ $\Delta t = 1$ s

Current actions

Scan HV Replay Zero Loudspeaker

Close Crystal calibration Help

Crystal calibration

0° (without target table)

Crystal: NaCl (2d = 564.0 pm)

Anode: Mo ($K_{\alpha} = 17.44$ keV)

Target angle: 7.2°

Progress

Target angle: 6.9° (-0.3°)
Sensor angle: 14.7° (+0.2°)
Counting rate: 2413 /s

The calibration of a crystal should be carried out after every crystal change or after every re-orientation of a crystal. This will be saved in the x-ray apparatus until the next calibration.

The 0° calibration (without target table) is sensible for Debye-Scherrer scans (powder targets).

For Bragg scans (single crystal), it is only necessary if the sensor angle has become so maladjusted that the automatic crystal calibration can no longer find the counting rate maximum.

Start search Cancel Adopt

Settings

X-Ray Apparatus Crystal General

Goniometer

Sensor

Target

Coupled $\beta =$

Parameters

U= 35.0 kV $\beta_{min} = 2.5^\circ$ $\Delta\beta = 0.1^\circ$

I= 1.00 mA $\beta_{max} = 30.0^\circ$ $\Delta t = 5$ s

Current actions

Scan HV Replay Zero Loudspeaker

Close Crystal calibration Help

35 kV, 30 kV, 25 kV

2.5°

30.0°


Passo Angular $\Delta\beta = 0.1^\circ$

Tempo de Contagem

1 mA

Bragg Plankc Transmission Moseley

$\beta / ^\circ$	$R_0 / 1/s$
15.3	56.4
15.4	56.0
15.5	56.6
15.6	52.4
15.7	52.8
15.8	51.2
15.9	51.0
16.0	45.2
16.1	48.6
16.2	40.2
16.3	45.0
16.4	49.0
16.5	44.2
16.6	39.2
16.7	44.4
16.8	40.4
16.9	38.4
17.0	38.8
17.1	37.6
17.2	43.4
17.3	33.0
17.4	34.0
17.5	32.8
17.6	29.6
17.7	33.0
17.8	32.6
17.9	31.6
18.0	37.4
18.1	35.0
18.2	25.4
18.3	32.0
18.4	33.2
18.5	26.8
18.6	27.0
18.7	28.6
18.8	27.4
18.9	27.4
19.0	29.0
19.1	26.6
19.2	30.0
19.3	30.0
19.4	38.2
19.5	47.0
19.6	47.2
19.7	40.4
19.8	28.2
19.9	24.6
20.0	24.8



Settings

X-Ray Apparatus Crystal General

Goniometer

Sensor

Target

Coupled $\beta = 20.0^\circ$

Parameter


U = 35.0 kV $\beta_{min} = 2.5^\circ$ $\Delta\beta = 0.1^\circ$


I = 1.00 mA $\beta_{max} = 20.0^\circ$ $\Delta t = 5$ s

Current actions:

Scan HV Replay Zero Loudspeaker


Close Crystal calibration Help





Bragg Planck Transmission Moseley

n_i / pm	R_0 / 1/s
24.6	2.8
25.6	2.8
26.6	2.8
27.6	3.0
28.5	3.2
29.5	1.8
30.5	2.8
31.5	3.6
32.5	5.6
33.4	4.4
34.4	4.0
35.4	6.0
36.4	4.0
37.4	6.6
38.4	10.2
39.3	24.4
40.3	55.2
41.3	85.8
42.3	131.4
43.3	168.4
44.3	205.4
45.2	247.8
46.2	278.0
47.2	324.8
48.2	369.4
49.2	391.6
50.1	410.6
51.1	449.0
52.1	457.2
53.1	465.4
54.1	502.8
55.0	495.2
56.0	529.6
57.0	514.8
58.0	546.0
59.0	531.6
59.9	529.4
60.9	608.2
61.9	798.6
62.9	1086.4
63.8	1100.0
64.8	779.4
65.8	574.4
66.8	533.8
67.8	529.4
68.7	734.4
69.7	1514.0
70.7	2196.2
71.7	2016.0



off

Settings

X-Ray Apparatus **Crystal** General

Lattice-plane spacing 2d: 564.0 pm (NaCl)

Representation

Energy conversion: off


Show all orders

Close Crystal calibration Help

X-Ray Apparatus

U = 35.0 kV
 I = 1.00 mA
 Δt = 5 s
 Δβ = 0.1°

0/s 40.0°
20.0°



Bragg Plank Transmission Moseley

$\beta / ^\circ$	$R_0 / 1/s$	R
15.3	56.4	1/s
15.4	56.0	

Properties

Delete Column

Select Character Size

Display Coordinates

Alt+C

Select Line Width

Show Values

 Show Connecting Lines

Select Rulers

 Show Grid

Logarithmic Representation

Zoom

Alt+Z

Zoom Off

Alt+O

Set Marker

Calculate Peak Center

Draw K-edge

Calculate Best-fit Straight Line

Calculate Straight Line through Origin

Calculate Integral

Delete Last Evaluation


Alt+Rück

Delete All Evaluations

Copy Table


Copy Diagram

Copy Window



Bragg Planck Transmission Moseley

n_i / pm	R_0 / 1/s
26.6	2.8
27.6	3.0
28.5	3.2
29.5	1.8
30.5	2.8
31.5	3.6
32.5	5.6
33.4	4.4
34.4	4.0
35.4	6.0
36.4	4.0
37.4	6.6
38.4	10.2
39.3	24.4
40.3	55.2
41.3	85.8
42.3	131.4
43.3	168.4
44.3	205.4
45.2	247.8
46.2	278.0
47.2	324.8
48.2	369.4
49.2	391.6
50.1	410.6
51.1	449.0
52.1	457.2
53.1	465.4
54.1	502.8
55.0	495.2
56.0	529.6
57.0	514.8
58.0	546.0
59.0	531.6
59.9	529.4
60.9	608.2
61.9	798.6
62.9	1086.4
63.8	1100.0
64.8	779.4
65.8	574.4
66.8	533.8
67.8	529.4
68.7	734.4
69.7	1514.0
70.7	2196.2
71.7	2016.8
72.6	1081.2
73.6	584.6



Settings

X-Ray Apparatus **Crystal** General

Lattice-plane spacing 2d: 564.0 pm (NaCl)

Representation

Energy conversion: off

Show all orders

Close Crystal calibration Help

X-Ray Apparatus

U = 35.0 kV
I = 1.00 mA
 $\Delta t = 5$ s
 $\Delta \beta = 0.1^\circ$

0/s 40.0°
20.0°

Properties

Delete Column

Select Character Size

Display Coordinates Alt+C

Select Line Width

Show Values

Show Connecting Lines

Select Rulers

Show Grid

Logarithmic Representation

Zoom Alt+Z

Zoom Off Alt+O

Set Marker

Calculate Peak Center

Draw K-edge

Calculate Best-fit Straight Line

Calculate Straight Line through Origin

Calculate Integral


Delete Last Evaluation Alt+Rück

Delete All Evaluations

Copy Table


Copy Diagram

Copy Window



Bragg Planck Transmission Moseley

n_i / pm	R_0 / 1/s
26.6	2.8
27.6	3.0
28.5	3.2
29.5	1.8
30.5	2.8
31.5	3.6
32.5	5.6
33.4	4.4
34.4	4.0
35.4	6.0
36.4	4.0
37.4	6.6
38.4	10.2
39.3	24.4
40.3	55.2
41.3	85.8
42.3	131.4
43.3	168.4
44.3	205.4
45.2	247.8
46.2	278.0
47.2	324.8
48.2	369.4
49.2	391.6
50.1	410.6
51.1	449.0
52.1	457.2
53.1	465.4
54.1	502.8
55.0	495.2
56.0	529.6
57.0	514.8
58.0	546.0
59.0	531.6
59.9	529.4
60.9	608.2
61.9	798.6
62.9	1086.4
63.8	1100.0
64.8	779.4
65.8	574.4
66.8	533.8
67.8	529.4
68.7	734.4
69.7	1514.0
70.7	2196.2
71.7	2016.8
72.6	1081.2
73.6	594.0



Valor do ângulo
 Incerteza do ângulo

Settings

X-Ray Apparatus Crystal General


Lattice-plane spacing 2d: [dropdown]

Representation

Energy conversion: off [dropdown]

Show all orders


Close Crystal calibration Help



$n_i = 70.93$ pm, $\sigma = 1.01$ pm

Bragg Plankc Transmission Moseley

$\beta / ^\circ$	$R_0 / 1/s$	R
2.5	2.8	1/s
2.6	2.8	
2.7	2.8	
2.8	3.0	
2.9	3.2	
3.0	1.8	
3.1	2.8	
3.2	3.6	
3.3	5.6	
3.4	4.4	
3.5	4.0	
3.6	6.0	
3.7	4.0	
3.8	6.6	
3.9	10.2	
4.0	24.4	
4.1	55.2	
4.2	85.8	
4.3	131.4	
4.4	168.4	
4.5	205.4	
4.6	247.8	
4.7	278.0	
4.8	324.8	
4.9	369.4	
5.0	391.6	
5.1	410.6	
5.2	449.0	
5.3	457.2	
5.4	465.4	
5.5	502.8	
5.6	495.2	
5.7	529.6	
5.8	514.8	
5.9	546.0	
6.0	531.6	
6.1	529.4	
6.2	608.2	
6.3	798.6	
6.4	1086.4	
6.5	1100.0	
6.6	779.4	
6.7	574.4	
6.8	533.8	
6.9	529.4	
7.0	734.4	
7.1	1514.0	
7.2	2196.2	



Repetir para U = 30 kV

Settings

X-Ray Apparatus Crystal General

Goniometer
 Sensor
 Target
 Coupled $\beta = 8.0^\circ$

Parameter
 U = 30.0 kV $\beta_{min} = 2.5^\circ$ $\Delta\beta = 0.1^\circ$
 I = 1.00 mA $\beta_{max} = 8.0^\circ$ $\Delta t = 3$ s

Current actions
 Scan HV Replay Zero Loudspeaker

Close Crystal calibration Help


X-Ray Apparatus

U = 30.0 kV
 I = 1.00 mA
 $\Delta t = 3$ s
 $\Delta\beta = 0.1^\circ$

1/s 16.0°
8.0°

Bragg Plankc Transmission Moseley

$\beta / ^\circ$	$R_0 / 1/s$	R
2.5	2.8	1/s
2.6	2.8	
2.7	2.8	
2.8	3.0	
2.9	3.2	
3.0	1.8	
3.1	2.8	
3.2	3.6	
3.3	5.6	
3.4	4.4	
3.5	4.0	
3.6	6.0	
3.7	4.0	
3.8	6.6	
3.9	10.2	
4.0	24.4	
4.1	55.2	
4.2	85.8	
4.3	131.4	
4.4	168.4	
4.5	205.4	
4.6	247.8	
4.7	278.0	
4.8	324.8	
4.9	369.4	
5.0	391.6	
5.1	410.6	
5.2	449.0	
5.3	457.2	
5.4	465.4	
5.5	502.8	
5.6	495.2	
5.7	529.6	
5.8	514.8	
5.9	546.0	
6.0	531.6	
6.1	529.4	
6.2	608.2	
6.3	798.6	
6.4	1086.4	
6.5	1100.0	
6.6	779.4	
6.7	574.4	
6.8	533.8	
6.9	529.4	
7.0	734.4	
7.1	1514.0	
7.2	2196.2	



Settings

X-Ray Apparatus Crystal General

Goniometer

Sensor

Target

Coupled $\beta = 8.0^\circ$

Parameter

U = 30.0 kV $\beta_{min} = 2.5^\circ$ $\Delta\beta = 0.1^\circ$

I = 1.00 mA $\beta_{max} = 8.0^\circ$ $\Delta t = 3$ s

Current actions

Scan HV Replay Zero Loudspeaker

Close Crystal calibration Help

X-Ray Apparatus


U = 30.0 kV
I = 1.00 mA
 $\Delta t = 3$ s
 $\Delta\beta = 0.1^\circ$

1/s 16.0°

8.0°

Bragg Plankc Transmission Moseley

$\beta / ^\circ$	$R_0 / 1/s$	R
2.5	2.8	1/s
2.6	2.8	
2.7	2.8	
2.8	3.0	
2.9	3.2	
3.0	1.8	
3.1	2.8	
3.2	3.6	
3.3	5.6	
3.4	4.4	
3.5	4.0	
3.6	6.0	
3.7	4.0	
3.8	6.6	
3.9	10.2	
4.0	24.4	
4.1	55.2	
4.2	85.8	
4.3	131.4	
4.4	168.4	
4.5	205.4	
4.6	247.8	
4.7	278.0	
4.8	324.8	
4.9	369.4	
5.0	391.6	
5.1	410.6	
5.2	449.0	
5.3	457.2	
5.4	465.4	
5.5	502.8	
5.6	495.2	
5.7	529.6	
5.8	514.8	
5.9	546.0	
6.0	531.6	
6.1	529.4	
6.2	608.2	
6.3	798.6	
6.4	1086.4	
6.5	1100.0	
6.6	779.4	
6.7	574.4	
6.8	533.8	
6.9	529.4	
7.0	734.4	
7.1	1514.0	
7.2	2196.2	



Repetir para U = 25 kV

Settings

X-Ray Apparatus Crystal General

Goniometer

Sensor

Target

Coupled $\beta = 8.0^\circ$

Parameter

U = 27.0 kV $\beta_{min} = 2.5^\circ$ $\Delta\beta = 0.1^\circ$

I = 1.00 mA $\beta_{max} = 8.0^\circ$ $\Delta t = 2$ s

Current actions

Scan HV Replay Zero Loudspeaker

Close Crystal calibration Help

X-Ray Apparatus

U = 27.0 kV

I = 1.00 mA

$\Delta t = 2$ s


$\Delta\beta = 0.1^\circ$

0/s 16.0°

8.0°

Bragg Plankc Transmission Moseley

$\beta / ^\circ$	$R_0 / 1/s$	R
2.5	2.8	1/s
2.6	2.8	
2.7	2.8	
2.8	3.0	
2.9	3.2	
3.0	1.8	
3.1	2.8	
3.2	3.6	
3.3	5.6	
3.4	4.4	
3.5	4.0	
3.6	6.0	
3.7	4.0	
3.8	6.6	
3.9	10.2	
4.0	24.4	
4.1	55.2	
4.2	85.8	
4.3	131.4	
4.4	168.4	
4.5	205.4	
4.6	247.8	
4.7	278.0	
4.8	324.8	
4.9	369.4	
5.0	391.6	
5.1	410.6	
5.2	449.0	
5.3	457.2	
5.4	465.4	
5.5	502.8	
5.6	495.2	
5.7	529.6	
5.8	514.8	
5.9	546.0	
6.0	531.6	
6.1	529.4	
6.2	608.2	
6.3	798.6	
6.4	1086.4	
6.5	1100.0	
6.6	779.4	
6.7	574.4	
6.8	533.8	
6.9	529.4	
7.0	734.4	
7.1	1514.0	
7.2	2195.2	



Importância de $\Delta\beta$ e Δt

Settings

X-Ray Apparatus Crystal General

Goniometer

Sensor

Target

Coupled $\beta = 2.5^\circ$

Parameter

U = 22.0 kV $\beta_{min} = 2.5^\circ$ $\Delta\beta = 0.2^\circ$

I = 1.00 mA $\beta_{max} = 8.0^\circ$ $\Delta t = 10$ s

Current actions

Scan HV Replay Zero Loudspeaker

Close Crystal calibration Help

X-Ray Apparatus

U = 22.0 kV

I = 1.00 mA

$\Delta t = 10$ s

$\Delta\beta = 0.2^\circ$

0/s

5.0°

2.5°



Bragg Plank Transmission Moseley

R₄ / 1/s R₅ / 1/s R


1/s

Selecione
a aba da coluna

Para deletar uma coluna
da tabela e o gráfico
correspondente



$R_4 / 1/s$	$R_5 / 1/s$	R
0.3	0.6	1/s
0.9	0.6	
0.6	0.6	
0.6	1.0	
0.6	0.4	
0.1	0.8	
0.6	0.0	
0.6	0.2	
0.8	1.2	
0.8	0.4	
0.8	0.8	
0.5	1.0	
0.8	0.2	
0.8	0.8	
0.4	0.6	
0.5	0.4	
0.5	0.8	
0.0	1.4	
0.6	0.4	
0.6	0.2	
0.7	1.0	
0.7	0.2	
0.7	0.4	
0.7	0.0	
0.5	0.2	
1.0	0.6	
1.0	0.6	
1.0	2.0	
0.6	4.2	
5.2	14.8	
34.5	30.2	
71.2	48.6	
85.2	65.2	
96.8	75.2	
102.6	94.0	
	98.4	
	109.4	



Settings

X-Ray Apparatus Crystal General

Goniometer
 Sensor
 Target
 Coupled $\beta = 8.0^\circ$

Parameter
 U = 22.0 kV $\beta_{min} = 2.5^\circ$ $\Delta\beta = 0.1^\circ$
 I = 1.00 mA $\beta_{max} = 8.0^\circ$ $\Delta t = 5$ s

Current actions
 Scan HV Replay Zero Loudspeaker

Close Crystal calibration Help


X-Ray Apparatus

Are you sure you want to delete table column 'R4'?

Sim Não


X-Ray Apparatus

U = 22.0 kV
 I = 1.00 mA
 $\Delta t = 5$ s
 $\Delta\beta = 0.1^\circ$



Bragg Planck Transmission Moseley


$\beta / ^\circ$	$R_0 / 1/s$	R
2.5	2.8	1/s
2.6	2.8	
2.7	2.8	
2.8	3.0	
2.9	3.2	
3.0	1.8	
3.1	2.8	
3.2	3.6	
3.3	5.6	
3.4	4.4	
3.5	4.0	
3.6	6.0	
3.7	4.0	
3.8	6.6	
3.9	10.2	
4.0	24.4	
4.1	55.2	
4.2	85.8	
4.3	131.4	
4.4	168.4	
4.5	205.4	
4.6	247.8	
4.7	278.0	
4.8	324.8	
4.9	369.4	
5.0	391.6	
5.1	410.6	
5.2	449.0	
5.3	457.2	
5.4	465.4	
5.5	502.8	
5.6	495.2	
5.7	529.6	
5.8	514.8	
5.9	546.0	
6.0	531.6	
6.1	529.4	
6.2	608.2	
6.3	798.6	
6.4	1086.4	
6.5	1100.0	
6.6	779.4	
6.7	574.4	
6.8	533.8	
6.9	529.4	
7.0	734.4	
7.1	1514.0	
7.2	2196.2	



Valores dos comprimentos de onda das transições K do Mo

	$\frac{E}{\text{keV}}$	$\frac{\nu}{\text{EHz}}$	$\frac{\lambda}{\text{pm}}$
K_{α}	17.443	4.2264	71.080
K_{β}	19.651	4.8287	63.095

$\text{keV} = 10^3 \text{ eV}$, $\text{EHz} = 10^{18} \text{ Hz}$, $\text{pm} = 10^{-12} \text{ m}$



A duração aproximada da tomada de dados do primeiro experimento é de 90 minutos (3 sessões de 30 minutos, uma para cada valor de U).

A duração aproximada da tomada de dados do segundo experimento é de 40 minutos (8 sessões de acordo com a tabela abaixo).

Antes de começar o 2º experimento, verifique se calculou todos os picos para os 3 valores de tensões do 1º experimento e para as três ordens de difração. Após isso, salve seus dados e apague todas as medidas para iniciar o 2º experimento com a tela limpa.



Bragg Planck Transmission Moseley

$\beta / ^\circ$ $R_0 / 1/s$ $R_1 / 1/s$ $R_2 / 1/s$ $R_3 / 1/s$ R

Properties

Delete Column

Select Character Size

Display Coordinates **Alt+C**

Select Line Width

Show Values

Show Connecting Lines

Select Rulers

Show Grid

Logarithmic Representation

Zoom Alt+Z

Zoom Off Alt+O

Set Marker

Calculate Peak Center

Draw K-edge

Calculate Best-fit Straight Line

Calculate Straight Line through Origin

Calculate Integral

Delete Last Evaluation **Alt+Rück**

Delete All Evaluations

Copy Table

Copy Diagram

Copy Window

7.1 1514.0 917.0 509.0 98.4

7.2 2196.2 1216.7 663.0 109.4

7.3 2016.0 1057.0 575.0 105.0

7.4 2016.0 1057.0 575.0 105.0

7.5 2016.0 1057.0 575.0 105.0

7.6 2016.0 1057.0 575.0 105.0

7.7 2016.0 1057.0 575.0 105.0

7.8 2016.0 1057.0 575.0 105.0

7.9 2016.0 1057.0 575.0 105.0

8.0 2016.0 1057.0 575.0 105.0

8.1 2016.0 1057.0 575.0 105.0

8.2 2016.0 1057.0 575.0 105.0

8.3 2016.0 1057.0 575.0 105.0

8.4 2016.0 1057.0 575.0 105.0

8.5 2016.0 1057.0 575.0 105.0

8.6 2016.0 1057.0 575.0 105.0

8.7 2016.0 1057.0 575.0 105.0

8.8 2016.0 1057.0 575.0 105.0

8.9 2016.0 1057.0 575.0 105.0

9.0 2016.0 1057.0 575.0 105.0

9.1 2016.0 1057.0 575.0 105.0

9.2 2016.0 1057.0 575.0 105.0

9.3 2016.0 1057.0 575.0 105.0

9.4 2016.0 1057.0 575.0 105.0

9.5 2016.0 1057.0 575.0 105.0

9.6 2016.0 1057.0 575.0 105.0

9.7 2016.0 1057.0 575.0 105.0

9.8 2016.0 1057.0 575.0 105.0

9.9 2016.0 1057.0 575.0 105.0

10.0 2016.0 1057.0 575.0 105.0

10.1 2016.0 1057.0 575.0 105.0

10.2 2016.0 1057.0 575.0 105.0

10.3 2016.0 1057.0 575.0 105.0

10.4 2016.0 1057.0 575.0 105.0

10.5 2016.0 1057.0 575.0 105.0

10.6 2016.0 1057.0 575.0 105.0

10.7 2016.0 1057.0 575.0 105.0

10.8 2016.0 1057.0 575.0 105.0

10.9 2016.0 1057.0 575.0 105.0

11.0 2016.0 1057.0 575.0 105.0

11.1 2016.0 1057.0 575.0 105.0

11.2 2016.0 1057.0 575.0 105.0

11.3 2016.0 1057.0 575.0 105.0

11.4 2016.0 1057.0 575.0 105.0

11.5 2016.0 1057.0 575.0 105.0

11.6 2016.0 1057.0 575.0 105.0

11.7 2016.0 1057.0 575.0 105.0

11.8 2016.0 1057.0 575.0 105.0

11.9 2016.0 1057.0 575.0 105.0

12.0 2016.0 1057.0 575.0 105.0

12.1 2016.0 1057.0 575.0 105.0

12.2 2016.0 1057.0 575.0 105.0

12.3 2016.0 1057.0 575.0 105.0

12.4 2016.0 1057.0 575.0 105.0

12.5 2016.0 1057.0 575.0 105.0

12.6 2016.0 1057.0 575.0 105.0

12.7 2016.0 1057.0 575.0 105.0

12.8 2016.0 1057.0 575.0 105.0

12.9 2016.0 1057.0 575.0 105.0

13.0 2016.0 1057.0 575.0 105.0

13.1 2016.0 1057.0 575.0 105.0

13.2 2016.0 1057.0 575.0 105.0

13.3 2016.0 1057.0 575.0 105.0

13.4 2016.0 1057.0 575.0 105.0

13.5 2016.0 1057.0 575.0 105.0

13.6 2016.0 1057.0 575.0 105.0

13.7 2016.0 1057.0 575.0 105.0

13.8 2016.0 1057.0 575.0 105.0

13.9 2016.0 1057.0 575.0 105.0

14.0 2016.0 1057.0 575.0 105.0

14.1 2016.0 1057.0 575.0 105.0

14.2 2016.0 1057.0 575.0 105.0

14.3 2016.0 1057.0 575.0 105.0

14.4 2016.0 1057.0 575.0 105.0

14.5 2016.0 1057.0 575.0 105.0

14.6 2016.0 1057.0 575.0 105.0

14.7 2016.0 1057.0 575.0 105.0

14.8 2016.0 1057.0 575.0 105.0

14.9 2016.0 1057.0 575.0 105.0

15.0 2016.0 1057.0 575.0 105.0

15.1 2016.0 1057.0 575.0 105.0

15.2 2016.0 1057.0 575.0 105.0

15.3 2016.0 1057.0 575.0 105.0

15.4 2016.0 1057.0 575.0 105.0

15.5 2016.0 1057.0 575.0 105.0

15.6 2016.0 1057.0 575.0 105.0

15.7 2016.0 1057.0 575.0 105.0

15.8 2016.0 1057.0 575.0 105.0

15.9 2016.0 1057.0 575.0 105.0

16.0 2016.0 1057.0 575.0 105.0

16.1 2016.0 1057.0 575.0 105.0

16.2 2016.0 1057.0 575.0 105.0

16.3 2016.0 1057.0 575.0 105.0

16.4 2016.0 1057.0 575.0 105.0

16.5 2016.0 1057.0 575.0 105.0

16.6 2016.0 1057.0 575.0 105.0

16.7 2016.0 1057.0 575.0 105.0

16.8 2016.0 1057.0 575.0 105.0

16.9 2016.0 1057.0 575.0 105.0

17.0 2016.0 1057.0 575.0 105.0

17.1 2016.0 1057.0 575.0 105.0

17.2 2016.0 1057.0 575.0 105.0

17.3 2016.0 1057.0 575.0 105.0

17.4 2016.0 1057.0 575.0 105.0

17.5 2016.0 1057.0 575.0 105.0

17.6 2016.0 1057.0 575.0 105.0

17.7 2016.0 1057.0 575.0 105.0

17.8 2016.0 1057.0 575.0 105.0

17.9 2016.0 1057.0 575.0 105.0

18.0 2016.0 1057.0 575.0 105.0

18.1 2016.0 1057.0 575.0 105.0

18.2 2016.0 1057.0 575.0 105.0

18.3 2016.0 1057.0 575.0 105.0

18.4 2016.0 1057.0 575.0 105.0

18.5 2016.0 1057.0 575.0 105.0

18.6 2016.0 1057.0 575.0 105.0

18.7 2016.0 1057.0 575.0 105.0

18.8 2016.0 1057.0 575.0 105.0

18.9 2016.0 1057.0 575.0 105.0

19.0 2016.0 1057.0 575.0 105.0

19.1 2016.0 1057.0 575.0 105.0

19.2 2016.0 1057.0 575.0 105.0

19.3 2016.0 1057.0 575.0 105.0

19.4 2016.0 1057.0 575.0 105.0

19.5 2016.0 1057.0 575.0 105.0

19.6 2016.0 1057.0 575.0 105.0

19.7 2016.0 1057.0 575.0 105.0

19.8 2016.0 1057.0 575.0 105.0

19.9 2016.0 1057.0 575.0 105.0

20.0 2016.0 1057.0 575.0 105.0

20.1 2016.0 1057.0 575.0 105.0

20.2 2016.0 1057.0 575.0 105.0

20.3 2016.0 1057.0 575.0 105.0

20.4 2016.0 1057.0 575.0 105.0

20.5 2016.0 1057.0 575.0 105.0

20.6 2016.0 1057.0 575.0 105.0

20.7 2016.0 1057.0 575.0 105.0

20.8 2016.0 1057.0 575.0 105.0

20.9 2016.0 1057.0 575.0 105.0

21.0 2016.0 1057.0 575.0 105.0

21.1 2016.0 1057.0 575.0 105.0

21.2 2016.0 1057.0 575.0 105.0

21.3 2016.0 1057.0 575.0 105.0

21.4 2016.0 1057.0 575.0 105.0

21.5 2016.0 1057.0 575.0 105.0

21.6 2016.0 1057.0 575.0 105.0

21.7 2016.0 1057.0 575.0 105.0

21.8 2016.0 1057.0 575.0 105.0

21.9 2016.0 1057.0 575.0 105.0

22.0 2016.0 1057.0 575.0 105.0

22.1 2016.0 1057.0 575.0 105.0

22.2 2016.0 1057.0 575.0 105.0

22.3 2016.0 1057.0 575.0 105.0

22.4 2016.0 1057.0 575.0 105.0

22.5 2016.0 1057.0 575.0 105.0

22.6 2016.0 1057.0 575.0 105.0

22.7 2016.0 1057.0 575.0 105.0

22.8 2016.0 1057.0 575.0 105.0

22.9 2016.0 1057.0 575.0 105.0

23.0 2016.0 1057.0 575.0 105.0


23.1 2016.0 1057.0 575.0 105.0

23.2 2016.0 1057.0 575.0 105.0

23.3 2016.0 1057.0 575.0 105.0


23.4 2016.0 1057.0 575.0 105.0

23.5 2016.0 1057.0 575.0 105.0



Bragg Planck Transmission Moseley

$\beta / ^\circ$	$R_0 / 1/s$	$R_1 / 1/s$	$R_2 / 1/s$	$R_3 / 1/s$
2.5	2.8	2.3	0.5	0.6
2.6	2.8	0.3	0.5	0.6
2.7	2.8	1.3	1.0	0.6
2.8	3.0	1.7	1.0	0.6
2.9	3.2	1.7	1.5	0.6
3.0	1.8	2.3	0.5	1.0
3.1	2.8	1.3	0.5	0.4
3.2	3.6	0.7	0.5	0.8
3.3	5.6	2.3	1.0	0.0
3.4	4.4	1.7	1.5	0.2
3.5	4.0	1.0	1.0	1.2
3.6	6.0	2.3	0.0	0.4
3.7	4.0	2.3	2.5	0.8
3.8	6.6	4.3	1.0	1.0
3.9	10.2	1.3	1.5	0.2
4.0	24.4	1.7	1.5	0.8
4.1	55.2	1.3	1.5	0.6
4.2	85.8	2.3	0.5	0.6
4.3	131.4	3.0	1.0	0.6
4.4	168.4	2.0	1.5	0.4
4.5	205.4	4.7	1.0	0.8
4.6	247.8	9.0	2.5	1.4
4.7	278.0	25.0	2.0	0.4
4.8	324.8	54.0	0.0	0.2
4.9	369.4	100.7	3.0	1.0
5.0	391.6	127.0	2.0	0.2
5.1	410.6	165.0	3.5	0.2
5.2	449.0	199.0	13.5	0.4
5.3	457.2	223.7	39.5	0.0
5.4	465.4	259.0	69.0	0.2
5.5	502.8	284.3	108.5	0.6
5.6	495.2	291.7	144.0	0.6
5.7	529.6	319.7	174.5	0.8
5.8	514.8	354.7	184.0	0.8
5.9	546.0	354.3	215.5	0.4
6.0	531.6	334.3	209.0	1.0
6.1	529.4	358.0	229.5	0.6
6.2	608.2	417.3	266.0	2.0
6.3	798.6	524.0	303.0	4.2
6.4	1086.4	627.0	360.5	14.8
6.5	1100.0	614.0	379.0	30.2
6.6	779.4	466.0	334.5	48.6
6.7	574.4	411.0	278.5	65.2
6.8	533.8	367.0	278.0	75.2
6.9	529.4	382.0	271.0	94.0
7.0	734.4	477.3	309.0	96.8
7.1	1514.0	917.0	509.0	98.4
7.2	2196.2	1216.7	663.0	109.4
7.3	2016.0	1057.0	575.0	105.0



Bragg Plank Transmission Moseley

 $\beta / ^\circ$ $R_0 / 1/s$ $R_1 / 1/s$ $R_2 / 1/s$ $R_3 / 1/s$ R 700

Properties

Delete Column

Select Character Size ▶

Display Coordinates

Alt+C

Select Line Width ▶

Show Values

✓ Show Connecting Lines

Select Rulers ▶

✓ Show Grid

Logarithmic Representation

Zoom

Alt+Z

Zoom Off

Alt+O

Set Marker ▶

Calculate Peak Center

Draw K-edge

Calculate Best-fit Straight Line

Calculate Straight Line through Origin

Calculate Integral

Delete Last Evaluation

Alt+Rück

Delete All Evaluations

Copy Table

Copy Diagram ▶

Copy Window

7.1	1514.0	917.0	509.0	98.4
7.2	2196.2	1216.7	663.0	109.4
7.3	2810.0	1057.0	575.0	105.0

2.5 3 3.5 4 4.5 5 5.5 6 6.5 7 7.5 8
 $\beta / ^\circ$

Bragg Plank Transmission Moseley

 $\beta / ^\circ$ $R_0 / 1/s$ $R_1 / 1/s$ $R_2 / 1/s$ $R_3 / 1/s$ R 700

Properties

Delete Column

Select Character Size

Display Coordinates

Alt+C

Select Line Width

Show Values

 Show Connecting Lines

Select Rulers

 Show Grid

Logarithmic Representation

Zoom

Alt+Z

Zoom Off

Alt+O

Set Marker

Calculate Peak Center

Draw K-edge

 Calculate Best-fit Straight Line

Calculate Straight Line through Origin

Calculate Integral

Delete Last Evaluation


Alt+Rück

Delete All Evaluations

Copy Table

Copy Diagram

Copy Window



7.1 1514.0 917.0 509.0 98.4

7.2 2196.2 1216.7 663.0 109.4


7.3 2010.0 1057.2 575.0 105.0

2.5 3 3.5 4 4.5 5 5.5 6 6.5 7 7.5 8 $\beta / ^\circ$ A = 372.0 / (s²), B = 3.9 °

© by LD Didactic GmbH, 1998-2008

Bragg Plankc Transmission Moseley

$\beta / ^\circ$	$R_0 / 1/s$	$R_1 / 1/s$	$R_2 / 1/s$	$R_3 / 1/s$
2.5	2.8	2.3	0.5	0.6
2.6	2.8	0.3	0.5	0.6
2.7	2.8	1.3	1.0	0.6
2.8	3.0	1.7	1.0	0.6
2.9	3.2	1.7	1.5	0.6
3.0	1.8	2.3	0.5	1.0
3.1	2.8	1.3	0.5	0.4
3.2	3.6	0.7	0.5	0.8
3.3	5.6	2.3	1.0	0.0
3.4	4.4	1.7	1.5	0.2
3.5	4.0	1.0	1.0	1.2
3.6	6.0	2.3	0.0	0.4
3.7	4.0	2.3	2.5	0.8
3.8	6.6	4.3	1.0	1.0
3.9	10.2	1.3	1.5	0.2
4.0	24.4	1.7	1.5	0.8
4.1	55.2	1.3	1.5	0.6
4.2	85.8	2.3	0.5	0.6
4.3	131.4	3.0	1.0	0.6
4.4	168.4	2.0	1.5	0.4
4.5	205.4	4.7	1.0	0.8
4.6	247.8	9.0	2.5	1.4
4.7	278.0	25.0	2.0	0.4
4.8	324.8	54.0	0.0	0.2
4.9	369.4	100.7	3.0	1.0
5.0	391.6	127.0	2.0	0.2
5.1	410.6	165.0	3.5	0.2
5.2	449.0	199.0	13.5	0.4
5.3	457.2	223.7	39.5	0.0
5.4	465.4	259.0	69.0	0.2
5.5	502.8	284.3	108.5	0.6
5.6	495.2	291.7	144.0	0.6
5.7	529.6	319.7	174.5	0.8
5.8	514.8	354.7	184.0	0.8
5.9	546.0	354.3	215.5	0.4
6.0	531.6	334.3	209.0	1.0
6.1	529.4	358.0	229.5	0.6
6.2	608.2	417.3	266.0	2.0
6.3	798.6	524.0	303.0	4.2
6.4	1086.4	627.0	360.5	14.8
6.5	1100.0	614.0	379.0	30.2
6.6	779.4	466.0	334.5	48.6
6.7	574.4	411.0	278.5	65.2
6.8	533.8	367.0	278.0	75.2
6.9	529.4	382.0	271.0	94.0
7.0	734.4	477.3	309.0	96.8
7.1	1514.0	917.0	509.0	98.4
7.2	2196.2	1216.7	663.0	109.4
7.3	2016.0	1057.2	575.0	105.2








Bragg Planck Transmission Moseley

U / kV λ_{min} / pm

Planck








Bragg Planck Transmission Moseley


n_i / pm	R_0 / 1/s	R_1 / 1/s	R_2 / 1/s	R_3 / 1/s
24.6	2.8	2.3	0.5	0.6
25.6	2.8	0.3	0.5	0.6
26.6	2.8	1.3	1.0	0.6
27.6	3.0	1.7	1.0	0.6
28.5	3.2	1.7	1.5	0.6
29.5	1.8	2.3	0.5	1.0
30.5	2.8	1.3	0.5	0.4
31.5	3.6	0.7	0.5	0.8
32.5	5.6	2.3	1.0	0.0
33.4	4.4	1.7	1.5	0.2
34.4	4.0	1.0	1.0	1.2
35.4	6.0	2.3	0.0	0.4
36.4	4.0	2.3	2.5	0.8
37.4	6.6	4.3	1.0	1.0
38.4	10.2	1.3	1.5	0.2
39.3	24.4	1.7	1.5	0.8
40.3	55.2	1.3	1.5	0.6
41.3	85.8	2.3	0.5	0.6
42.3	131.4	3.0	1.0	0.6
43.3	168.4	2.0	1.5	0.4
44.3	205.4	4.7	1.0	0.8
45.2	247.8	9.0	2.5	1.4
46.2	278.0	25.0	2.0	0.4
47.2	324.8	54.0	0.0	0.2
48.2	369.4	100.7	3.0	1.0
49.2	391.6	127.0	2.0	0.2
50.1	410.6	165.0	3.5	0.2
51.1	449.0	199.0	13.5	0.4
52.1	457.2	223.7	39.5	0.0
53.1	465.4	259.0	69.0	0.2
54.1	502.8	284.3	108.5	0.6
55.0	495.2	291.7	144.0	0.6
56.0	529.6	319.7	174.5	0.8
57.0	514.8	354.7	184.0	0.8
58.0	546.0	354.3	215.5	0.4
59.0	531.6	334.3	209.0	1.0
59.9	529.4	358.0	229.5	0.6
60.9	608.2	417.3	266.0	2.0
61.9	798.6	524.0	303.0	4.2
62.9	1086.4	627.0	360.5	14.8
63.8	1100.0	614.0	379.0	30.2
64.8	779.4	466.0	334.5	48.6
65.8	574.4	411.0	278.5	65.2
66.8	533.8	367.0	278.0	75.2
67.8	529.4	382.0	271.0	94.0
68.7	734.4	477.3	309.0	96.8
69.7	1514.0	917.0	509.0	98.4
70.7	2196.2	1216.7	663.0	109.4
71.7	2016.0	1057.0	575.0	105.0






Bragg Plank Transmission Moseley

n_i / pm	R_0 / 1/s	R_1 / 1/s	R_2 / 1/s	R_3 / 1/s
24.6	2.8	2.3	0.5	0.6
25.6	2.8	0.3	0.5	0.6
26.6	2.8	1.3	1.0	0.6
27.6	3.0	1.7	1.0	0.6
28.5	3.2	1.7	1.5	0.6
29.5	1.8	2.3	0.5	1.0
30.5	2.8	1.3	0.5	0.4
31.5	3.6	0.7	0.5	0.8
32.5	5.6	2.3	1.0	0.0
33.4	4.4	1.7	1.5	0.2
34.4	4.0	1.0	1.0	1.2
35.4	6.0	2.3	0.0	0.4
36.4	4.0	2.3	2.5	0.8
37.4	6.6	4.3	1.0	1.0
38.4	10.2	1.3	1.5	0.2
39.3	24.4	1.7	1.5	0.8
40.3	55.2	1.3	1.5	0.6
41.3	85.8	2.3	0.5	0.6
42.3	131.4	3.0	1.0	0.6
43.3	168.4	2.0	1.5	0.4
44.3	205.4	4.7	1.0	0.8
45.2	247.8	9.0	2.5	1.4
46.2	278.0	25.0	2.0	0.4
47.2	324.8	54.0	0.0	0.2
48.2	369.4	100.7	3.0	1.0
49.2	391.6	127.0	2.0	0.2
50.1	410.6	165.0	3.5	0.2
51.1	449.0	199.0	13.5	0.4
52.1	457.2	223.7	39.5	0.0
53.1	465.4	259.0	69.0	0.2
54.1	502.8	284.3	108.5	0.6
55.0	495.2	291.7	144.0	0.6
56.0	529.6	319.7	174.5	0.8
57.0	514.8	354.7	184.0	0.8
58.0	546.0	354.3	215.5	0.4
59.0	531.6	334.3	209.0	1.0
59.9	529.4	358.0	229.5	0.6
60.9	608.2	417.3	266.0	2.0
61.9	798.6	524.0	303.0	4.2
62.9	1086.4	627.0	360.5	14.8
63.8	1100.0	614.0	379.0	30.2
64.8	779.4	466.0	334.5	48.6
65.8	574.4	411.0	278.5	65.2
66.8	533.8	367.0	278.0	75.2
67.8	529.4	382.0	271.0	94.0
68.7	734.4	477.3	309.0	96.8
69.7	1514.0	917.0	509.0	98.4
70.7	2196.2	1216.7	663.0	109.4
71.7	2016.0	1057.0	575.0	105.0





Bragg Plank Transmission Moseley

U / kV	λ_{\min} / pm
35.0	38.9
30.0	45.4
27.0	50.9
22.0	61.8

