

The principle and method of ELISA

Principle and method of the experiment

- [Antibody](#)
- [Qualitative and quantitative measurements of proteins using antibodies](#)
- [Fractionation and purification of proteins](#)

Featured

- [Campaign](#)
- [New products](#)
- [Antibody search](#)
- [Product list](#)

Contact


Enzyme-linked immunosorbent assay (ELISA) is a method of target antigen (or antibody) capture in samples using a specific antibody (or antigen), and of target molecule detection/quantitation using an enzyme reaction with its substrate.

The principle

In ELISA, various antigen-antibody combinations are used, always including an enzyme-labeled antigen or antibody, and enzyme activity is measured colorimetrically.



The enzyme activity is measured using a substrate that changes color when modified by the enzyme. Light absorption of the product formed after substrate addition is measured and converted to numeric values. Depending on the antigen-antibody combination, the assay is called a direct ELISA, indirect ELISA, sandwich ELISA, competitive ELISA etc.

ELISA: An example of an assay using a 96-well plate.





The yellow color indicates that the target protein is present. The higher degree of the color, the higher concentration of the target protein.

Direct ELISA



A target protein (or a target antibody) is immobilized on the surface of microplate wells and incubated with an enzyme-labeled antibody to the target protein (or a specific antigen to the target antibody). After washing, the activity of the microplate well-bound enzyme is measured.


Indirect ELISA



A target protein is immobilized on the surface of microplate wells and incubated with an antibody to the target protein (the primary antibody), followed by a secondary antibody against the primary antibody. After washing, the activity of the microplate well-bound enzyme is measured.

Although indirect ELISA requires more steps than direct ELISA, labeled secondary antibodies are commercially available, eliminating the need to label the primary antibody.

Sandwich ELISA




An antibody to a target protein is immobilized on the surface of microplate wells and incubated first with the target protein and then with another target protein-specific antibody, which is labeled with an enzyme. After washing, the activity of the microplate well-bound enzyme is measured. The immobilized antibody (orange) and the enzyme-labeled antibody (green) must recognize different epitopes of the target protein.

★AZ's Point★



Compared to direct ELISA, the sandwich ELISA (combining antibodies to two different epitopes on the target protein) has a higher specificity. Sandwich ELISA is useful for applications that require a high accuracy.

Competitive ELISA



An antibody specific for a target protein is immobilized on the surface of microplate wells and incubated with samples containing the target protein and a known amount of enzyme-labeled target protein. After the reaction, the activity of the microplate well-bound enzyme is measured.

When the antigen level in the sample is high, the level of antibody-bound enzyme-labeled antigen is lower and the color is lighter. Conversely, when it is low, the level of antibody-bound enzyme-labeled antigen is higher and the color, darker. The graph above and to the right illustrates the correlation between absorption and antigen levels in samples.

★AZ's Point★



When a target antigen is a small molecule, such as histamine, pesticide, and dioxin, two antibodies cannot simultaneously bind to the antigen in sandwich ELISA. Competitive ELISA is useful for the measurement of low molecular weight targets.

Immobilization of antibody/antigen on microplate wells

How are proteins bound (adsorbed) to the wells?

Immobilization is commonly mediated by hydrophobic interaction or covalent bonds. Microplates that are pre-activated for different immobilization methods are commercially available. The choice of method of immobilization depends on the type of protein to be immobilized. In general, covalent bonding is preferred for smaller molecules, such as peptides, and hydrophobic interaction (physical adsorption) is preferred for larger molecules, such as proteins. An appropriate method should be determined for each assay system.




Pre-treatment of microplates

Many types of ELISA microplates are commercially available, including the hydrophobic type, hydrophilic type as well as types surface-activated with amino or carboxyl groups. Surface-activated microplates are useful for covalent bonding.


It is important to select the correct type of ELISA microplate for the intended use.

ELISA procedure

Sandwich ELISA ※An example performed at MBL


Step-by-step procedure	
<p>Preparation of reagents and equipment</p>	 <p>96-well plate ▼ ▼ Reservoir</p> <p>▲ 12-Channel pipette</p>
<p>Immobilization of antibody Add diluted antibody to each well of a 96-well ELISA plate. Seal the plate to prevent evaporation, and allow it to incubate at 4°C for 15-18 hours to immobilize the antibody.</p>	
<p>Washing Remove the diluted antibody, and wash 3 times with washing solution.</p>	
<p>Blocking Add blocking buffer to each well, and allow it to incubate at 37°C for 1 hour to reduce non-specific binding of the target protein to the well.</p>	
<p>Washing Remove the blocking buffer, and wash 3 times with washing solution.</p>	
<p>Addition of samples Dilute the samples with sample dilution buffer, and add 100 µL of each sample to each well. For the calibration curve, prepare a dilution series of the standard on the same plate. Allow it to incubate at 37°C for 1 hour.</p>	
<p>Washing Remove the samples, and wash 5 times with washing solution.</p>	
<p>Addition of the detection antibody</p>	

Dilute the detection antibody in sample dilution buffer, and add 100 μ L to each well.
Allow it to incubate at 37°C for 1 hour.





Washing
After reaction, remove the detection antibody, and wash 5 times with washing solution.

Addition of enzyme-labeled secondary antibody
Dilute an enzyme-labeled secondary antibody with sample dilution buffer, and add 100 μ L to each well.
Allow it to incubate at 37°C for 1 hour.



Washing
After reaction, remove the secondary antibody, and wash 5 times with washing solution.


Add a substrate solution.
Allow it to incubate as the color develops.




Dilution series of a standard for the calibration curve

Add a stop solution to stop the reaction when the color is sufficiently developed.

Adding stop solution ▶




▼ After the addition of stop solution





Dilution series of a standard for the calibration curve

Measure the absorption at 450 nm by a plate reader.

Plate reader ▼



[Related topics] Links to ELISA-related products

Antibodies suitable for ELISA	List of primary antibodies List of secondary antibodies
ELISA kits	Kits to make sandwich ELISA: Ab-Match Kits ELISA kits and related products for drug discovery  ELISA kits and related products for metabolism research 
When Blocking	Bovine Serum Albumin Reagent Grade Powder (pH7.0)

[Next page: The principle and method of immunoprecipitation \(IP\)](#)

[Previous page: The principle and method of Western blotting \(WB\)](#)

Related links

Antibody

Antibody basics

- [What are antibodies?](#)
- [Structure of antibodies](#)

[The role of antibodies](#)[Types of antibodies](#)

Antibodies as a research tool

[How to generate antibodies](#)[How to select antibodies](#)[Labeled antibodies](#)[How to label antibodies](#)[Main causes of non-specific reactions](#)[How to reduce non-specific reactions](#)[Tags and Tag Antibodies](#)

Qualitative and quantitative measurements of proteins using antibodies

[The principle and method of Western blotting \(WB\)](#)[The principle and method of ELISA](#)[The principle and method of immunoprecipitation \(IP\)](#)[The principle and method of co-immunoprecipitation \(Co-IP\)](#)

Fractionation and purification of proteins

[The principle and method of polyacrylamide gel electrophoresis \(SDS-PAGE\)](#)[The principle and method of chromatography](#)[HOME](#)[Security & Privacy](#) | [About this site](#)[↑ TOP](#)

Product category

<Research area>

- Autophagy
- Immunology
- Stem cell research
- RNA-RNP network
- Epigenetics
- Cancer
- Metabolism
- Drug discovery
- Apoptosis
- Ubiquitin proteasome
- Neuroscience
- Multiple transmembrane protein
- Circadian rhythm
- Cell surface antigens
- Cytokine & Growth Factors
- Signal transduction
- Cellular stress
- Exosomes

<Tools>

- Epitope tags
- MHC tetramer
- Fluorescent proteins
- Control antibodies
- Immunoglobulin
- Cell culture reagents
- Flow cytometry reagents
- 2nd antibodies
- Latex particles
- Magnetic particles

<Brands>

- CycLex
- CircuLex
- Amalgaam
- CoralHue™
- HRP-Direct
- RiboCluster Profiler™
- IMMUTEX™
- Magnosphere™

Custom services

information

- New products
- Campaign
- Discontinued Products
- Topics
- Exhibitions
- Literatures offer

Learn

- Principle and method of the experiment
- Webinar

Download

- Consent and Request forms
- Brochures & leaflets
- MBL RUO LETTER
- Immunology News
- Pathway posters

Contact

MBL MEDICAL & BIOLOGICAL LABORATORIES CO., LTD.
A JSR Life Sciences Company

Copyright © 2017 MEDICAL & BIOLOGICAL LABORATORIES CO., LTD. All rights reserved.