

Clonogenic Assays Protocols

Clonogenic assays serve as a useful tool to test whether a given cancer therapy can reduce the clonogenic survival of tumor cells. A colony is defined as a cluster of at least 50 cells that can often only be determined microscopically. A clonogenic assay is the method of choice to determine cell reproductive death after treatment with ionizing radiation, but can also be used to determine the effectiveness of other cytotoxic agents.

Crystal Violet Assay

Materials and Reagents

- Cell culture medium (DMEM 10% FCS)
- Phosphate buffered saline (PBS)
- Trypsin/ EDTA
- Crystal violet
- Fixation solution
- Cell culture petri dishes or six-well plates
- Hemocytometer
- Stereomicroscope (e.g., Nikon Eclipse, model: TS100)
- CO₂ Incubator

Cell preparation:

1. Culture CF41 and MeLn cells according to their requirements.
2. Remove medium, and then rinse cells with 10 ml PBS.
3. Add 4 ml 0.25% trypsin to the cells and incubate for 1-5 min until the cells appear round and detach from the substrate.
4. Add 10 ml medium with 10% FBS, and detach the cells by pipetting.
5. Count the cells using a hemocytometer.
6. Note: It is critical to get a relatively accurate number for the cells.
7. Prepare desired seeding concentration: 1×10^3 to 4×10^4 cells per milliliter of complete culture medium, and then seed cell into 6mm dishes.

Assay setup:

1. Incubate cells for a few hours to overnight in a CO₂ incubator at 37 °C and allow them to attach to the plate/dish.
2. Treat the cells as necessary with chemicals, radiation or a combination of both.
3. Incubate the cells in a CO₂ incubator at 37 °C for 1-2 weeks until cells in control plates have formed colonies that are of a substantially good size (50 cells per colony is the minimum for scoring).

Fixation and staining:

1. Remove medium, and then rinse cells with 10 ml PBS.
2. Remove PBS and add 2-3 ml of fixation solution and leave the dishes/plates at room temperature (RT) for 5 min.
3. Remove fixation solution.
4. Add 0.5% crystal violet solution and incubate at RT for 2 h.
5. Remove crystal violet carefully and wash the dishes/plates in distilled water to rinse off crystal violet.
6. Air-dry the dishes/plates on a table cloth at RT for up to a few days.

Data analysis:

1. Count number of colonies with a stereomicroscope.
2. Calculate plating efficiency (PE) and surviving fraction (SF).
3. $PE = \text{no. of colonies formed} / \text{no. of cells seeded} \times 100\%$
4. $SF = \text{no. of colonies formed after treatment} / \text{no. of cells seeded} \times PE$

Soft agar assay (Clonogenic anchorage-independent assay)

1. Prepare single-cell suspensions according to the Cell Preparation item up to step 6;
2. Mix 1×10^3 to 4×10^4 cells in 1 mL of 0.6% agarose prepared in complete medium;
3. Plate cell-agarose mixes in culture dishes coated with 5mm of 1% agarose (already prepared by your instructors);
4. Incubate plates at 37 °C and 5% CO₂ for 2 weeks to allow colony formation;
5. The number of colonies formed and the colony area (mm²) are determined under the light microscope;
6. Analyze using ImageJ software (NIH, Bethesda, MD, USA).